

INSTITUTE OF MEDICAL BIOCHEMISTRY AND LABORATORY MEDICINE

Introduction to Work in Laboratory

Measuring volumes, filtration, centrifugation,
solubility, separation

Practical in Medical Biochemistry
General Medicine

Martin Vejražka



2022/23

Content

TASK 1 – SEPARATION OF DISPERSION BY MEANS OF CENTRIFUGATION	3
1. PRINCIPLE	3
2. PROCEDURE	4
3. TASKS	4
TASK 2 – SEPARATION OF A DISPERSION BY SHAKING WITH CHLOROFORM	4
1. PRINCIPLE	4
2. PROCEDURE	5
3. TASKS	5
TASK 3 – ADSORPTION OF DISSOLVED COMPOUND ON ACTIVATED CHARCOAL	5
1. PRINCIPLE	5
2. PROCEDURE	7
3. TASKS	7
TASK 4 – SUBLIMATION OF CAFFEINE FROM COFFEE	7
1. PRINCIPLE	7
2. PROCEDURE	8
3. TASKS	8
TASK 5 – DISSOLUTION OF A WEAK ACID IN AN ENVIRONMENT OF VARIOUS PH	9
1. PRINCIPLE	9
2. PROCEDURE	9
3. TASKS	9

Task 1 – Separation of dispersion by means of centrifugation

1. Principle

In this practical class, diluted milk coloured with two dyes will be used as a model dispersion. The first dye is reddish lipophilic substance, named Sudan III (Fig. 1); the second one is a hydrophilic bluish food-dye, brilliant blue (also known under the code E-133, Fig. 2). Milk is an emulsion of the type oil in water. The hydrophilic component forms the dispersion environment and the dispersed part is created mostly by fat. This emulsion is stable especially thanks to proteins that possess both the hydrophilic and lipophilic parts.

In the first experiment, we will try to separate components of emulsion by means of simple centrifugation. Two different centrifugation forces will be used: a bench-top microcentrifuge reaches approximately 300×g, a high-speed centrifuge will be set to 30,000×g. The time of centrifugation will be the same in both cases.




Fig. 1: Sudan III




Fig. 2: Brilliant blue (E-133)

2. Procedure

1. Milk
2. Purified water
3. Sudan III 4 g/l in alcohol
4. Brilliant blue (E-133) 0.25 g/l in purified water

Procedure

Sample preparation

Dilute milk with purified water in the proportion 1:2 – in a beaker, mix 5 mL of milk with 10 mL of purified water. Use the glass pipette to measure volumes. Mix gently.

Mix approx. 10 mL of diluted milk, 10 drops of red colorant Sudan III in alcohol, and 10 drops of aqueous solution of the blue colorant brilliant blue (E-133). Seal the test-tube with laboratory film (Parafilm) and shake with vortex.

Centrifugation

1. Measure 1 mL of the mixture to each of two micro-test tubes. Close the tubes with caps.
2. One test-tube is to be spun with a bench-top microcentrifuge for 10 minutes.
3. The second test-tube is to be centrifuged for 10 minutes at 30,000×g in a high-speed centrifuge.

3. Tasks

Compare separation of the mixture and explain the differences.

Task 2 – Separation of a dispersion by shaking with chloroform

1. Principle

A mixture of compounds of various polarities can be split between two immiscible liquid phases. Chloroform is a strongly hydrophobic, non-polar organic solvent. It is heavier than water. Chloroform is practically immiscible with hydrophilic substances. When mixtures of hydrophilic and hydrophobic compounds (emulsion of fat in water in our case) are shaken with chloroform, hydrophobic components are dissolved in chloroform while the hydrophilic ones remain in the aqueous phase.




Fig. 3: Separation of emulsion by shaking with chloroform.

Left: emulsion stained with blue hydrophilic and red hydrophobic dye during shaking with chloroform. Right: after phase separation.

2. Procedure

1. Diluted milk stained with Sudan III and brilliant blue from the previous experiment
2. Chloroform (trichlormethan) 

Procedure

1. Pour approx. 4 mL of stained diluted milk from the previous experiment into an equilibration test-tube. Add approximately the same volume of chloroform.
2. Close the test-tube and shake with vortex (the highest speed) three times, 30 seconds each time. The content must always be completely shuffled.
3. Put the closed test-tube in a rack. Let it stay until the aqueous (upper) phase is separated from the organic (bottom) phase.

Caution: Residues that contain chloroform are disposed of as hazardous waste.

3. Tasks

1. Compare the appearance of two phases. Compare the colour with the original mixture. Explain the differences.
2. With the structural formula of the two dyes, explain which parts of their molecules are responsible for the observed behaviour.

Task 3 – Adsorption of dissolved compound on activated charcoal

1. Principle

Activated charcoal is pure carbon in a form that has an extremely large surface area. It consists of graphite particles, in which micropores (about the size of the order of several nanometres to several tens of nanometres). Porosity of particles greatly increases their adsorptive surface.

Activated charcoal binds, adsorbs, many compounds in a non-specific way. It is generally reported that the adsorption of substances on the carbon increases with molecular weight, polarity of molecules, the number of multiple bonds and functional groups.

Activated charcoal is often used for purification of solutions from small quantities of contaminants. It is usually added to the solution as a powder, the resulting suspension is stirred over a period of time in order to bind contaminants, and carbon particles are then filtered or centrifuged out. Both steps are combined in this experiment as we will filter the solution through a membrane containing activated charcoal.



Fig. 4: Filter assembly with a charcoal filter.




Fig. 5: Adsorption on activated charcoal

Right: a residue of the mixture after separation by shaking with chloroform from the previous experiment. Left: the bluish dye is adsorbed to activated charcoal during filtration through a charcoal filter. Pink color of the filtrate is caused by traces of the red hydrophobic dye that contaminates the hydrophilic phase and that is much less adsorbed to activated charcoal.

2. Procedure

1. Separated mixture from the previous experiment
2. Charcoal filters (Schleicher&Schuller No. 508, diameter 9 cm), equipment for filtration

Procedure

1. Pre-wet one circle of charcoal filter with purified water. Fold it and put it to a funnel.
Caution: Activated charcoal stains.
2. Aspire the aqueous phase from the previous experiment with a dropper and filter it through the charcoal filter. The filtrate is to be collected in a clean test-tube.

3. Tasks

Compare the appearance of the filtrate with the original mixture and with the aqueous phase from the previous experiment. Compare it with the diluted milk as well. Explain the differences.

Task 4 – Sublimation of caffeine from coffee

1. Principle

Some compounds can be separated from a mixture by means of distillation or sublimation. The fact that the purified compound has the boiling point different from other compounds in the mixture is used in both cases. When the mixture is heated the isolated compound is converted into gaseous state. The gas is derived from the mixture and cooled down so that it condensates back to a liquid or solid state. In case the compound is isolated from a liquid, the process is described as distillation. In sublimation, the solid state is converted directly to gas when the mixture is heated.

In this experiment we will isolate caffeine from roasted and ground coffee. Caffeine is relatively thermo-stable and its boiling point is 178 °C. Its melting point is higher (about 235 °C), therefore sublimation takes place easily when it is heated. When the vapour is cooled down, white acicular crystals of caffeine are formed.




Fig 6: Caffeine




Fig. 7: Sublimation of caffeine from coffee

2. Procedure

1. Roasted and ground coffee desiccated at 120 °C
2. Sea sand, aluminium foil, beaker 250 mL, watch glass

Procedure

1. Desiccate approx. 7 g of coffee for 5 minutes at 120 °C.
2. Prepare a sand bath: Pour the sea sand into a bowl of thicker aluminum foil. Place on the cooker in the fume hood and heat to about 200 ° C (stage 3).
3. Pour the desiccated ground coffee to a clean and thoroughly polished beaker. Put the beaker into the sand and cover it with a polished watch glass.
4. Let the sublimation run for about 5 minutes. Then take the beaker with glass from the sand and let it cool down.

Caution: Both the beaker and the watch glass are hot! Use forceps or an insulating glove.

3. Tasks

Observe crystals of caffeine on the wall of the beaker.




Fig. 8: Crystals of caffeine

Task 5 – Dissolution of a weak acid in an environment of various pH

1. Principle

Polarity of organic acids and bases strongly depends on pH of environment. Carboxylic acids, for example, form carboxylate anions in alkaline solutions and thus turn strongly polar; they become therefore easily soluble in aqueous solution. On the other hand, in acidic environment, carboxylic acids get protonated, less polar and their solubility is lower. This is important for example for absorption of drugs: most drugs behave as weak acids or weak bases, and pH differs in various sections of gastrointestinal tract.

In this experiment solubility of salicylic acid in acidic, basic and approximately neutral environment will be examined.






Figure 9: Salicylic acid

2. Procedure

1. Crystalline salicylic acid
2. Hydrochloric acid 2 mol/l 
3. Sodium hydroxide 2 mol/l 

Procedure

1. Pour a small amount (about 0.1 g) of crystalline salicylic acid to each of three test-tubes.
2. Add approx. 2 mL of purified water to the first test-tube, 2 mL of diluted hydrochloric acid to the second one, and 2 mL of diluted sodium hydroxide to the third one.
3. Mix the content of test-tubes and observe how quickly crystals of salicylic acid dissolve.
4. When crystals are completely dissolved in the test-tube with sodium hydroxide, add approx. 3 mL of diluted hydrochloric acid and mix carefully. Watch changes in the mixture.

Caution: The mixture may get hot during neutralisation. It may start to boil and sprinkle. Work so that the mouth of the test-tube is facing away of you and people standing around!

3. Tasks

Compare solubility of salicylic acid in acidic, neutral and basic environment. Explain what happened after neutralisation of the basic solution of the compound. Use chemical equations to describe the observed changes.