

Instructions for the practical lesson on biochemistry

Topic: Calcium, phosphorus, sodium, metabolism of bone tissue

Task 1: Estimation of total calcium in serum and urine

Reagents:

Commercial kit *Bio-La-Test CALCIUM Liquid 250*, made by Erba-Lachema s.r.o.



1. **Working solution:** MES buffer pH 6.5 100 mmol/l, arsenazo III 200 μmol/l
2. **Standard solution of calcium 2.5 mmol/l**
3. **Serum** – unknown sample (infectious material)
4. **Urine** – unknown sample (infectious material) **diluted 4×**

Procedure:

N.B.: Ultrapure deionized water must be used for this experiment both for the blank and for rinsing the cuvettes in between measurements.

Measure reagents to **plastic (!) single-use test tubes** according to the following table:

Measure in ml:	Serum sample (Tube No. 1)	Urine sample (Tube No. 2)	Standard (Tube No. 3)	Blank (Tube No. 4)
Working solution	1.0	1.0	1.0	1.0
Serum	0.01	-	-	-
Urine	-	0.01	-	-
Standard solution	-	-	0.01	-
Ultrapure water	-	-	-	0.01

Mix and incubate 1 minute at 37 °C. Measure the absorbances of serum, urine and standard against the blank at 650 nm in new 1 cm plastic cuvette. The coloration is stable for one hour.

Calculations:

Concentration of total calcium in the serum (S-Ca):

$$\text{S-Calcium (mmol/l)} = \frac{A_{\text{serum}}}{A_{\text{standard}}} \times C_{\text{standard}}$$

Concentration of calcium in the urine (U-Calcium):

$$\text{U-Calcium (mmol/l)} = \frac{A_{\text{urine}}}{A_{\text{standard}}} \times C_{\text{standard}} \times \text{Dilution of urine}$$

Daily output of calcium into urine (dU-Calcium):

$$\text{dU-Calcium (mmol/24 hrs)} = \text{U-Calcium (mmol/l)} \times \text{Volume of urine (liters/24 hrs)}$$

Task 2: Estimation of inorganic phosphate in serum and urine

Reagents:

Commercial kit Bio-La-Test PHOSPHORUS UV Liquid 250, made by Erba-Lachema s.r.o.



1. **Working solution:** Sulfuric acid 210 mmol/l, ammonium molybdate 650 mmol/l
1. **Standard solution of phosphate 1.61 mmol/l**
2. **Serum** – unknown sample (infectious material)
3. **Urine** – unknown sample (infectious material) **diluted 25×**

Procedure:

N.B.: Ultrapure deionized water must be used for this experiment both for the blank and for rinsing the cuvettes in between measurements.

Measure reagents to **plastic (!) single-use test tubes** according to the following table:

Measure in ml:	Serum sample (Tube No. 1)	Urine sample (Tube No. 2)	Standard (Tube No. 3)	Blank (Tube No. 4)
Working solution	1.0	1.0	1.0	1.0
Serum	0.01	-	-	-
Urine (diluted)	-	0.01	-	-
Standard solution	-	-	0.01	-
Ultrapure water	-	-	-	0.01

Mix and incubate 5 minutes at 37 °C. Measure the absorbances of serum, urine and standard against the blank at 340 nm in new 1 cm plastic cuvette. The coloration is stable for one hour.

Calculations:

Concentration of inorganic phosphate in serum (fS-P inorg.):

$$\text{fS-Inorganic phosphate (mmol/l)} = \frac{A_{\text{serum}}}{A_{\text{standard}}} \times C_{\text{standard}}$$

Concentration of inorganic phosphate in urine (U-P inorg.):

$$\text{U- Inorganic phosphate (mmol/l)} = \frac{A_{\text{urine}}}{A_{\text{standard}}} \times C_{\text{standard}} \times \text{Dilution of urine}$$

Daily output of inorganic phosphate into urine (dU-P):

$$\text{dU-Inorg. phosphate (mmol/24 hrs)} = \text{U-Inorg. phosphate (mmol/l)} \times \text{Volume of urine (liters/24 hrs)}$$

Task 3: Estimation of catalytic concentration of alkaline phosphatase and its isoforms

Reagents:

Commercial kit *Bio-La-Test Alkalická fosfatáza 120* manufactured by Erba-Lachema s.r.o.. is used.

- Buffer:** N-methyl-D-glucamine pH 10.1, 0.35 mol/l, NaCl 70.0 mmol/l, MgCl₂ 0.5 mmol/l
- Substrate:** 4-nitrophenyl phosphate disodium salt 15 mmol/l in the buffer described above
(Concentrations in the final reaction mixture)



- Inhibitor:** Sodium hydroxide 1 mol/l, Chelaton 3 30 mmol/l
- Serum** – unknown sample (infectious material)

Procedure:

Two serum aliquots are prepared for this experiment: 'S' (native) and 'T' (subjected to heat inactivation at 56 °C for 10 minutes).

Proceed according to the table:

Measure in ml:	Native serum sample (S1)	Control for native sample (S2)	Heat-inactivated sample (T1)	Control for heat-inactivated (T2)
Buffer	1.0	1.0	1.0	1.0
Native serum (S)	0.02	-	-	-
Heat-inactivated serum (T)	-	-	0.02	-
Mix and incubate 5 minutes at 37 °C				
Substrate	0.2	0.2	0.2	0.2
Mix and incubate exactly 10 minutes at 37 °C				
Inhibitor	0.5	0.5	0.5	0.5
Native serum (S)	-	0.02	-	-
Heat-inactivated serum (T)	-	-	-	0.02
Mix and within 30 minutes measure the absorbances of S1, S2, T1, and T2 against deionized water at 420 nm in 1 cm cuvette.				

Calculations:

Subtraction of control absorbances:

$$\Delta A_S = A_{S1} - A_{S2} \qquad \Delta A_T = A_{T1} - A_{T2}$$

Total catalytic concentration of ALP:

Calculation is based on the molar absorption coefficient of 4-nitrophenol.

$$\text{Total ALP } (\mu\text{kat/l}) = \Delta A_S \times 10.263$$

Catalytic concentration of the liver isoenzyme:

Calculation assumes that 100 % of the bone and one third of the liver isoenzyme have been inactivated by heat.



$$\text{Liver isoenzyme ALP } (\mu\text{kat/l}) = 1.5 \times \Delta A_T \times 10.263$$

Catalytic concentration of the bone isoenzyme:

$$\text{Bone isoenzyme ALP } (\mu\text{kat/l}) = \text{Total ALP } (\mu\text{kat/l}) - \text{Liver isoenzyme ALP } (\mu\text{kat/l})$$

Task 4 Solubility of various calcium salts

Reagents:

1. **Calcium chloride**, CaCl_2 , solid substance
2. **Calcium carbonate**, CaCO_3 , solid substance
3. **Calcium phosphate**, $\text{Ca}_2(\text{PO}_4)_3$, solid substance
4. **Ammonium oxalate**, $(\text{NH}_4)_2(\text{COO})_2$, saturated solution 
5. **Hydrochloric acid**, HCl , diluted 
6. **Sodium hydrogen carbonate**, NaHCO_3 , 10 g/l
7. **Na₂EDTA**, saturated solution
8. **Lactose**, 10 g/l

A. Solubility of calcium salts in water and HCl

Procedure:

- a. Take four short glass test tubes and pour a pinch of solid calcium chloride (tube No. 1), calcium carbonate (tube No. 2) or calcium phosphate (tubes No. 3 and 4).
- b. Add about 1.5 ml of deionized water to each tube.
- c. Mix and observe whether the different calcium salts dissolve in water.
- d. To the tubes No. 2, 3, and 4 try to add diluted hydrochloric acid drop wise until the precipitate just dissolves; record the number of drops to the table in your report.
- e. Next, to the tube No. 3 containing dissolved $\text{Ca}_3(\text{PO}_4)_2$ add sodium hydrogen carbonate solution drop wise and observe formation of a precipitate.
- f. To the tube No. 4 add about 8 drops of saturated Na_2EDTA solution followed by the sodium hydrogen carbonate solution (the same number of drops as to the tube No. 3). Na_2EDTA prevents precipitation.
- g. Summarize all your observations to the table in your report.

B. Influence of some food components on solubility of calcium salts

Procedure:

Prepare two mixtures according to the table and observe whether a precipitate of calcium salt results.

	Tube No. 1	Tube No. 2
CaCl ₂	1 measure	1 measure
Deionized water	1 ml	1 ml
Ammonium oxalate	cca 5 drops	-
Lactose	-	cca 5 drops

Task 5 Quantitative estimation of Na⁺ in urine

Reagents:

1. **Dilution buffer:** 1.98 g of acetic acid 100 % and 4.09 g of ethanolamine in 1,000 ml, pH alkaline.
2. **Standard NaCl solutions** 0.1 mol/l, 0.01 mol/l, and 0.001 mol/l in the dilution buffer.
3. **Urine sample** of unknown Na⁺ concentration, **diluted 10×** with the dilution buffer.

Procedure:

- a. Dilute the urine sample 10× with the buffer: mix 5 ml of sample with 45 ml of buffer.
- b. Follow the provided instructions for use of pH-meter as a mV-meter and measure the electrode potentials in the calibration solutions of NaCl 1 – 3 as well as in the diluted urine sample. Do not rinse the electrode with distilled water in between measurements; wipe it gently (avoid rubbing!) instead. Record all measured values of E (mV) into table in your report.
- c. Create a calibration graph from the measured values of the electrode potentials and the corresponding pNa values of the standard solutions. Plot the pNa values on the x axis and the potential in mV on the y axis. Next, use the calibration graph to read the pNa for the analyzed urine sample.
- d. Convert the pNa value of diluted urine sample into Na⁺ concentration:

$$\text{pNa} = -\log[\text{Na}^+] \quad [\text{Na}^+]_{\text{diluted urine}} = 10^{-\text{pNa}}$$

- e. Correct for dilution of urine:

$$\text{U-Na}^+ = 10 \times [\text{Na}^+]_{\text{diluted urine}}$$

Further calculations:

Daily output of Na⁺ into urine:

$$dU\text{-Na}^+ = U\text{-Na}^+ (\text{mol/l}) \times \text{Vol. urine (liters/24 hrs)}$$

Fractional excretion (FE) of Na⁺:

$$FE_{Na} = \frac{U_{Na} \times P_{Cr}}{U_{Cr} \times P_{Na}}$$

U_{Cr} Creatinine in urine (mmol/l)

P_{Cr} Creatinine in serum (mmol/l)

U_{Na} Sodium in urine (mmol/l)

P_{Na} Sodium in serum (mmol/l)

Tubular resorption (TR) of Na⁺:

$$TR_{Na} = 1 - FE_{Na}$$