

Date Name Group

Instructions and lab report form for the practical lesson on biochemistry

Topic: Separation techniques

Task 1: Separation of hemoglobin and potassium ferricyanide using gel filtration

Objectives:

- Measure absorption spectrum of ferricyanide to identify wavelength that can be used for its spectrophotometric determination.
- Construct calibration curve for ferricyanide.
- Explain why you cannot measure concentration of both analytes in unknown sample directly, using simple spectrophotometry? (*Check the absorption spectrum of both compounds.*)
- Separate both substances using gel filtration and determine their concentration.
- Determine elution volume of both compounds.

Principle:

Observation/Results:

Calibration:

Ferricyanide, $\lambda =$

$c_1 =$	$A =$
$c_2 =$	$A =$
$c_3 =$	$A =$
$c_4 =$	$A =$
$c_5 =$	$A =$
$c_6 =$	$A =$

Conclusion:

Task 2: Thin-layer chromatography of plant pigments


Objectives:

- Separate a mixture of plant pigments using thin-layer chromatography with two different mobile phases (highly non-polar hexane and a partially polar hexane: acetone mixture).
- Identify individual pigments on the chromatogram, developed in the partially polar mobile phase (your teacher will help you). Calculate R_f values for all pigments.
- Identify pigments on the chromatogram, developed in hexane (by yourself, without assistance from your teacher).
- Decide, which of the green spots belongs to chlorophyll a and which to chlorophyll b. Explain your conclusion. Refer to structural formulas of chlorophylls below.
- Highlight the difference(s) between lutein and β-carotene (by referring to the structural formulas below) and decide which one is more polar. Does it explain their different mobility in the partially polar mobile phase?
- Explain why most of the pigments do not separate in the non-polar mobile phase.


Principle:

Observation:


Conclusion:




Chlorophyll a



Chlorophyll b



β-carotene



Lutein

Task 3: Dialysis

Objectives: Demonstrate semi-permeability of dialysis membrane.

Principle:

Observation:

Conclusion: