

Date Name Group

Lab report from the practical lesson on biochemistry

Topic: Metabolism of calcium and phosphate. Osteoporosis

Task 1: Estimation of inorganic phosphate in serum and urine

Principle:

Results:

	Serum sample (Tube No. 1)	Urine sample (Tube No. 2)	Standard (Tube No. 3)	Blank (Tube No. 4)
A 340 nm				-

Calculations:

Concentration of inorganic phosphate in serum (fS-P inorg.):

$$\text{fS-Inorganic phosphate (mmol/l)} = \frac{A_{\text{serum}}}{A_{\text{standard}}} \times C_{\text{standard}}$$

$$\text{fS-Inorganic phosphate (mmol/l)} = \frac{\dots\dots\dots}{\dots\dots\dots} \times \dots\dots\dots = \dots\dots\dots$$

Concentration of inorganic phosphate in urine (U-P inorg.):

$$\text{U-Inorganic phosphate (mmol/l)} = \frac{A_{\text{urine}}}{A_{\text{standard}}} \times C_{\text{standard}} \times \text{Dilution of urine}$$

$$\text{U-Inorganic phosphate (mmol/l)} = \frac{\dots\dots\dots}{\dots\dots\dots} \times \dots\dots\dots \times \dots\dots\dots = \dots\dots\dots$$

Daily output of inorganic phosphate into urine (dU-P):

$$\text{dU-Inorg. phosphate (mmol/24 hrs)} = \text{U-Inorg. phosphate (mmol/l)} \times \text{Volume of urine (liters/24 hrs)}$$

$$\text{dU-Inorg. phosphate (mmol/24 hrs)} = \dots\dots\dots \times \dots\dots\dots = \dots\dots\dots$$

Conclusion:

Evaluate both serum phosphate and daily output of phosphate to urine.

Task 2: Estimation of total calcium in serum and urine**Principle:****Results:**

	Serum sample (Tube No. 1)	Urine sample (Tube No. 2)	Standard (Tube No. 3)	Blank (Tube No. 4)
A 650 nm				-

Calculations:***Concentration of total calcium in the serum (S-Ca):***

$$\text{S-Calcium (mmol/l)} = \frac{A_{\text{serum}}}{A_{\text{standard}}} \times C_{\text{standard}}$$

$$\text{S-Calcium (mmol/l)} = \frac{\dots\dots\dots}{\dots\dots\dots} \times \dots\dots\dots = \dots\dots\dots$$

Concentration of calcium in the urine (U-Calcium):

$$\text{U- Calcium (mmol/l)} = \frac{A_{\text{urine}}}{A_{\text{standard}}} \times C_{\text{standard}} \times \text{Dilution of urine}$$

$$\text{U- Calcium (mmol/l)} = \frac{\dots\dots\dots}{\dots\dots\dots} \times \dots\dots\dots \times \dots\dots\dots = \dots\dots\dots$$

Daily output of calcium into urine (dU-Calcium):

$$\text{dU-Calcium (mmol/24 hrs)} = \text{U-Calcium (mmol/l)} \times \text{Volume of urine (liters/24 hrs)}$$

$$\text{dU-Calcium (mmol/24 hrs)} = \dots \times \dots = \dots$$

Conclusion:

Evaluate both serum calcium and daily output of calcium to urine.

Task 3: Estimation of catalytic concentration of alkaline phosphatase and its isoforms**Principle:**

What is chemistry of the ALP estimation? How is separate measurement of different isoenzymes achieved?

A. End-point assay with estimation of isoforms**Results and calculations:**

	Native serum sample (S1)	Control for native sample (S2)	Heat-inactivated sample (T1)	Control for heat-inactivated (T2)
A 420 nm				

Subtraction of control absorbances:

$$\Delta A_S = A_{S1} - A_{S2}$$

$$\Delta A_S = \dots - \dots = \dots$$

$$\Delta A_T = A_{T1} - A_{T2}$$

$$\Delta A_T = \dots - \dots = \dots$$

Total catalytic concentration of ALP:

$$\text{Total ALP } (\mu\text{kat/l}) = \Delta A_S \times 10.263 = \dots \times 10.263 = \dots$$

Catalytic concentration of the liver isoenzyme:

$$\text{Liver isoenzyme ALP } (\mu\text{kat/l}) = 1.5 \times \Delta A_T \times 10.263 = \dots$$

Catalytic concentration of the bone isoenzyme:

$$\text{Bone isoenzyme ALP } (\mu\text{kat/l}) = \text{Total ALP } (\mu\text{kat/l}) - \text{Liver isoenzyme ALP } (\mu\text{kat/l})$$

$$\text{Bone isoenzyme ALP } (\mu\text{kat/l}) = \dots - \dots = \dots$$

B. Total ALP – kinetic assay

Results and calculations:

Measured absorbances:	Difference per minute:
A ₀	
A ₁	$\Delta A_1 = A_1 - A_0 = \dots$
A ₂	$\Delta A_2 = A_2 - A_1 = \dots$
A ₃	$\Delta A_3 = A_3 - A_2 = \dots$
A ₄	$\Delta A_4 = A_4 - A_3 = \dots$
A ₅	$\Delta A_5 = A_5 - A_4 = \dots$

Mean difference of absorbance per minute $\Delta A_{420}/\text{min.}$:

$$\text{Catalytic concentration ALP } (\mu\text{kat/l}) = \Delta A_{420}/\text{min.} \times 72.8 = \dots$$

Conclusion:

Evaluate both total ALP and the proportion of both isoforms.

Task 4: Solubility of various calcium salts

Principle:

A. Solubility of calcium salts in water and HCl

Results:

	Tube No. 1 CaCl ₂	Tube No. 2 CaCO ₃	Tube No. 3 Ca ₃ (PO ₄) ₂	Tube No. 4 Ca ₃ (PO ₄) ₂ + Na ₂ EDTA
Solubility in water				
Solubility in HCl	–			
Solubility in NaHCO ₃	–	–		

B. Influence of some food components on solubility of calcium salts

Results:

	Tube No. 1 CaCl ₂ + Ammonium oxalate	Tube No. 2 CaCl ₂ + Lactose
Result		

Conclusion: